

CE  
2292

**SURGICLEAN®**

**Absorbable Gauze of  
Oxidized  
Regenerated  
Cellulose**

*Control Bleeding Improve Healing*

Hangzhou Singclean Medical Products Co., Ltd.  
Add: No. 125(E), 10th Street, Hangzhou Qiantang  
New Area, Zhejiang, China 310018  
P.C.: 310018  
Tel: +86-571-63431868 Fax: +86-571-63431886  
E-mail: [marketing@hzhc.com](mailto:marketing@hzhc.com)  
Web: [www.singclean.net](http://www.singclean.net)

More information on



## What is Surgiclean® Absorbable Gauze of Oxidized Regenerated Cellulose?

It's medical gauze for surgery.

The raw material is regenerated cellulose made from specialty fiber through oxidation and trimming. This smart hemostatic material can be decomposed and absorbed during the healing process of human wound tissue.

As a supporting hemostatic device, Surgiclean® can be used for various surgical procedures. It is especially ideal in capillary, vein and arteries where the traditional hemostasis is limited such as ligation. At the same time, it can be used instead of ordinary gauze during surgery to stop bleeding as well as reduce the possibility of postoperative adhesions.

### ● How it works

#### Step1


Blood pass through the gauze slowly

#### Step2

Gauze absorb the blood and forms into a tight physical seal at the bleeding surface

#### Step3

Hemostasis within 3 minutes, fully absorbed with 7-30 days.



### storage

In clean, cool, dry and ventilated environment

Brand	Specification
<b>SURGICLEAN®</b>	2 × 3cm
	3 × 5cm
	5 × 5cm
	5 × 7.5cm
	5 × 8cm
	5 × 10cm
	5 × 25cm
	5 × 30cm
	5 × 35cm

It can be customizable



## Application

Applicable Departments:

General Surgery / Nail Surgery / Hepatobiliary And Pancreatic Surgery / Gastrointestinal Surgery / Gynecology

Various Open And Endoscopic Surgery To Stop Bleeding

Applicable Surgery Type: Appendectomy / Breast Surgery / Hernia Repair / Gastrectomy / Hepatectomy / Hepatobiliary Bladder Surgery Gynecological Surgery, Thyroid Surgery, Treatment Of Superficial Injuries, Etc.

**Ready to use, right out of the box**

## WHY? SURGICLEAN®

Ordinary look with extraoRdinary performanCe

Plant-based material

Low PH creates bactericidal environment

Multiple sizes/shapes

Sterile and Biodegradable

Hemostasis within 3 minutes

Absorbed with 7-30 days